

NM interference in the DCF assay

Quenching effects – DCF

AUTHORED BY:	DATE:
Cordula Hirsch	20.01.14

REVIEWED BY:	DATE:
Harald Krug	09.04.2014

APPROVED BY:	DATE:

DOCUMENT HISTORY

Effective Date	Date Revision Required	Supersedes
15.02.2014	DD/MM/YYYY	DD/MM/YYYY

Version	Date	Description of the Change		Author / Changed by
1.0	15/02/2014	All	Initial Document	Cordula Hirsch
1.1	20/10/2016	1	NM "solvent" for pre-dilutions (old: ddH ₂ O; new: HBSS; paragraph 6.2)	Cordula Hirsch

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		1/11

Table of Content

1	Introduction.....	3
2	Principle of the Method	3
3	Applicability and Limitations	3
4	Related Documents	3
5	Equipment and Reagents	4
5.1	Equipment	4
5.2	Reagents	4
5.3	Reagent Preparation	4
5.3.1	1x concentrated Hank’s Balanced Salt Solution (HBSS)	4
5.3.2	2x concentrated HBSS	5
5.3.3	DCF.....	5
5.3.4	Pluronic F-127.....	5
6	Procedure	6
6.1	Flow chart	6
6.2	Dilution of nanomaterials.....	6
6.3	Preparation of DCF working solution:	8
6.4	Distribution into 96-well plate	8
6.5	Application of NMs and measurement	9
6.6	Data evaluation	10
7	Quality Control, Quality Assurance, Acceptance Criteria.....	10
8	Health and Safety Warnings, Cautions and Waste Treatment.....	10
9	Abbreviations	10
10	References.....	11

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		2/11

1 Introduction

NMs have been shown to interfere with several colorimetric as well as fluorimetric read outs (e.g. Belyanskaya, 2007; Casey, 2007; Guo, 2008; Monteiro-Riviere, 2006; Pulskamp, 2007; Wörle-Knirsch, 2006). Especially carbon-based NMs may quench existing fluorescence signals (for a review see Kroll et al., 2009). The detection of potential quenching effects of NMs in the DCF-assay is addressed in this SOP.

2 Principle of the Method

The 2', 7'-dichlorodihydrofluorescein diacetate (H₂DCF-DA) assay is a widely used *in vitro* ROS-detection method. The non-fluorescent dye (H₂DCF-DA) is a chemically reduced form of fluorescein and cell-permeable. Intracellular esterases cleave off the diacetate (DA) moiety which renders the molecule (H₂DCF) sensitive to oxidation by ROS. In its oxidized form dichlorofluorescein (DCF) is highly fluorescent and easily detectable e.g. using a fluorescent plate reader.

The fluorescent DCF molecule is commercially available. It is used in a cell free 96-well plate setup to uncover potential quenching effects of NMs.

3 Applicability and Limitations

The results from these cell free controls cannot be calculated against values from cellular measurements. They serve as qualitative estimations of NM only reactions that do not involve cellular contribution.

4 Related Documents

Table 1: Documents needed to proceed according to this SOP and additional NM-related interference control protocols.

Document ID	Document Title
O_DCF_A549	<i>Detection of reactive oxygen species in A549 cell – DCF assay in A549 cells</i>
O_DCF_THP-1	<i>Detection of reactive oxygen species in THP-1 cell – DCF assay in THP-1 cells</i>
M_NM suspension_metal oxides	<i>Suspending and diluting Nanomaterials – Metal oxides and NM purchased as monodisperse suspensions</i>
M_NM suspension_carbon based	<i>Suspending and diluting Nanomaterials – Carbon based nanomaterials</i>

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		3/11

5 Equipment and Reagents

5.1 Equipment

- Flat bottom 96-well cell culture plates
- 15 ml conical tubes (polypropylene or polystyrene; e.g. from Falcon)
- Vortex®
- Multichannel pipette (with 12 positions; volume range per pipetting step at least from 50 µl to 200 µl)
- Fluorescence reader for multi-well plates (to measure excitation/emission at wavelength maxima of: $\lambda_{ex}=485$ nm and $\lambda_{em}=528$ nm)

5.2 Reagents

- 2',7'-Dichlorofluorescein (DCF) [CAS number: 76-54-0]
- Calcium chloride dihydrate ($\text{CaCl}_2 \cdot 2 \text{H}_2\text{O}$) [CAS number: 10035-04-8]
- D-Glucose [CAS number: 50-99-7]
- Disodium phosphate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$) [CAS number: 10028-24-7]
- Ethanol [CAS number: 64-17-5]
- Magnesium chloride hydrate ($\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$) [CAS number: 7791-18-9]
- Magnesium sulfate heptahydrate ($\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$) [CAS number: 7487-88-9]
- Pluronic F-127 [CAS number: 9003-11-6]
- Pluronic F-127 [CAS number: 9003-11-6]
- Potassium chloride (KCl) [CAS number: 7447-40-7]
- Potassium hydrogen phosphate (KH_2PO_4) [CAS number: 7778-77-0]
- Sodium chloride (NaCl) [CAS number: 8028-77-1]
- Sodium hydrogen carbonate (NaHCO_3) [CAS number: 7542-12-3]

5.3 Reagent Preparation

5.3.1 1x concentrated Hank's Balanced Salt Solution (HBSS)

1 g/l	D-glucose
185 mg/l	$\text{CaCl}_2 \cdot 2 \text{H}_2\text{O}$
400 mg/l	KCl
60 mg/l	KH_2PO_4
100 mg/l	$\text{MgCl}_2 \cdot 6 \text{H}_2\text{O}$
100 mg/l	$\text{MgSO}_4 \cdot 7 \text{H}_2\text{O}$
8 g/l	NaCl
350 mg/l	NaHCO_3
60 mg/l	$\text{Na}_2\text{HPO}_4 \cdot 2 \text{H}_2\text{O}$

Dissolve all reagents in ddH₂O and adjust the pH to 7.4. Store at 4°C.

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		4/11

5.3.2 2x concentrated HBSS

2 g/l	D-glucose
370 mg/l	CaCl ₂ * 2 H ₂ O
800 mg/l	KCl
120 mg/l	KH ₂ PO ₄
200 mg/l	MgCl ₂ * 6 H ₂ O
200 mg/l	MgSO ₄ * 7 H ₂ O
16 g/l	NaCl
700 mg/l	NaHCO ₃
120 mg/l	Na ₂ HPO ₄ * 2 H ₂ O

Dissolve all reagents in ddH₂O and adjust the pH to 7.4. Store at 4°C.

5.3.3 DCF

Stock:

- 500 µM in Ethanol: 0.2 mg/ml
- Prepare aliquots and store at -20°C. **Note:** can be stored for several month.

Working concentration:

- 2.5 µM in HBSS: dilute 15 µl of the 500 µM stock solution in 3 ml HBSS

5.3.4 Pluronic F-127

Stock:

- 160 ppm in ddH₂O: 160 µg/ml (=16 mg/100 ml)

Can be stored for up to 4 weeks at 4°C.

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		5/11

6 Procedure

6.1 Flow chart




Figure 1: Brief outline of the workflow.

6.2 Dilution of nanomaterials

For this SOP we distinguish two types of nanomaterials (NM) according to their solvent, suspension properties and highest concentrations used in the assay. See also respective related documents (3).

- (1) Metal oxide NM, Polystyrene beads and all NM delivered as monodisperse suspensions by the supplier: solvent either determined by the supplier or ddH₂O; sub-diluted in ddH₂O; highest concentration in assay 100 µg/ml
- (2) Carbon based NM: suspended and sub-diluted in 160 ppm Pluronic F-127; highest concentration in assay 80 µg/ml

Volumes given in the following dilution schemes are enough for one 96-well plate.

Note: “Mixing” in the context of diluting NMs means, the solvent containing tube is put on a continuously shaking Vortex® and the previous sub-dilution (or stock suspension, respectively) is put dropwise into the shaking solvent. The resulting suspension stays on the Vortex® for additional 3 seconds before proceeding with the next sub-dilution.

Note: A series of experiments showed that quenching reactions of certain NM (i.e. gold nanoparticles) can be stronger in ddH₂O compared to HBSS as a solvent. As cell based experiments are run in HBSS, this interference protocol was adapted accordingly to mimic the cellular situation as good as possible. Thus all NM pre-dilutions described below use HBSS as a solvent.

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		6/11

(1) Metal oxide NM:

Prepare serial sub-dilutions of the stock suspension (1 mg/ml) in HBSS.

- Label nine conical tubes (15 ml total volume) with 1 to 9 (relates to steps 1-9 below).
 - Mix 1.5 ml 2x concentrated HBSS and 300 µl ddH₂O in tube no. 1 (mixture A).
 - Add 1.5 ml (1x concentrated) HBSS to tubes 2 to 9.
1. 1.2 ml of 1 mg/ml stock suspension are mixed with 1.8 ml of mixture A → 400 µg/ml (1)
 2. 1.5 ml of 400 µg/ml (1) are mixed with 1.5 ml HBSS → 200 µg/ml (2)
 3. 1.5 ml of 200 µg/ml (2) are mixed with 1.5 ml HBSS → 100 µg/ml (3)
 4. 1.5 ml of 100 µg/ml (3) are mixed with 1.5 ml HBSS → 50 µg/ml (4)
 5. 1.5 ml of 50 µg/ml (4) are mixed with 1.5 ml HBSS → 25 µg/ml (5)
 6. 1.5 ml of 25 µg/ml (5) are mixed with 1.5 ml HBSS → 12.5 µg/ml (6)
 7. 1.5 ml of 12.5 µg/ml (6) are mixed with 1.5 ml HBSS → 6.25 µg/ml (7)
 8. 1.5 ml of 6.25 µg/ml (7) are mixed with 1.5 ml HBSS → 3.13 µg/ml (8)
 9. 1.5 ml of 3.13 µg/ml (8) are mixed with 1.5 ml HBSS → 1.56 µg/ml (9)

As these sub-dilutions will be mixed 1:2 with the DCF dye (or ddH₂O) in the 96-well plate (see Figure 2 and Figure 3), the final concentrations will be halved to the following final concentrations:

1. 200 µg/ml
2. 100 µg/ml
3. 50 µg/ml
4. 25 µg/ml
5. 12.5 µg/ml
6. 6.25 µg/ml
7. 3.13 µg/ml
8. 1.56 µg/ml
9. 0.78 µg/ml
10. Solvent control: HBSS

(2) Carbon based NM:

Prepare serial sub-dilutions of the stock suspension (500 µg/ml) in 160 ppm Pluronic F-127:

- Label nine conical tubes (15 ml total volume) with 1 to 9 (relates to steps 1-9 below).
- Add 1.5 ml 2x concentrated HBSS to tube no. 1.
- Mix 10 ml 2x HBSS with 10 ml 160 ppm Pluronic F-127.
Note: Using this mixture (A) for NM sub-dilutions will result in 1x HBSS containing the appropriate amount of NM as well as Pluronic F-127.
- Add 1.08 ml mixture (A) to tube no. 2.
- Add 1.5 ml mixture (A) to tubes 3 to 9.

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		7/11

1. 1.5 ml of 500 µg/ml stock suspension are mixed with 1.5 ml 2x HBSS → 250 µg/ml (1)
2. 1.92 ml of 250 µg/ml (1) are mixed with 1.08 ml (A) → 160 µg/ml (2)
3. 1.5 ml of 160 µg/ml (2) are mixed with 1.5 ml (A) → 80 µg/ml (3)
4. 1.5 ml of 80 µg/ml (3) are mixed with 1.5 ml (A) → 40 µg/ml (4)
5. 1.5 ml of 40 µg/ml (4) are mixed with 1.5 ml (A) → 20 µg/ml (5)
6. 1.5 ml of 20 µg/ml (5) are mixed with 1.5 ml (A) → 10 µg/ml (6)
7. 1.5 ml of 10 µg/ml (6) are mixed with 1.5 ml (A) → 5 µg/ml (7)
8. 1.5 ml of 5 µg/ml (7) are mixed with 1.5 ml (A) → 2.5 µg/ml (8)
9. 1.5 ml of 2.5 µg/ml (8) are mixed with 1.5 ml (A) → 1.25 µg/ml (9)

As these sub-dilutions will be mixed 1:2 with the DCF dye (or ddH₂O) in the 96-well plate (see Figure 2 and Figure 3), the final concentrations will be halved to the following final concentrations:

1. 125 µg/ml
2. 80 µg/ml
3. 40 µg/ml
4. 20 µg/ml
5. 10 µg/ml
6. 5 µg/ml
7. 2.5 µg/ml
8. 1.25 µg/ml
9. 0.625 µg/ml
10. Solvent control: mixture (A) = 80 ppm Pluronic F-127

6.3 Preparation of DCF working solution:

Dilute the 500 µM stock solution 1:200 in ddH₂O resulting in a 2.5 µM working solution. For one 96-well plate a final volume of 3 ml is needed.

- Mix 2985 µl ddH₂O with 15 µl of the 500 µM stock solution.

6.4 Distribution into 96-well plate

- Add 50 µl ddH₂O into each of the outermost wells (black wells in Figure 2, A1-A12; H1-H12; B1-G1; B12-G12).
- Add 50 µl ddH₂O also into wells B2-D11 (also black in Figure 2).
- Add 50 µl 2.5 µM DCF working solution into each blue well (Figure 2, E2 to G11).

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		8/11




Figure 2: Plate layout.

Rows A to D as well as outermost wells (E1-H1; E12-H12; H2-H11) receive 50 μ l ddH₂O each. Inner wells (E2-G11) receive 50 μ l 2.5 μ M DCF.

6.5 Application of NMs and measurement

Note: All NM dilutions have to be vortexed directly before application.

- Add NM sub-dilutions (see 6.2 “Dilution of nanomaterials”) to appropriate wells according to Figure 3 and the table below:

wells	metal oxide NM concentration	carbon based NM concentration
B11-G11	400 μ g/ml	250 μ g/ml
B10-G10	200 μ g/ml	160 μ g/ml
B9-G9	100 μ g/ml	80 μ g/ml
B8-G8	50 μ g/ml	40 μ g/ml
B7-G7	25 μ g/ml	20 μ g/ml
B6-G6	12.5 μ g/ml	10 μ g/ml
B5-G5	6.25 μ g/ml	5 μ g/ml
B4-G4	1.13 μ g/ml	2.5 μ g/ml
B3-G3	1.56 μ g/ml	1.25 μ g/ml
B2-D2	solvent (HBSS)	solvent (80 ppm Pluronic F-127 in HBSS)




Figure 3: Distribution of NM dilutions. 50 μ l of the NM sub-dilutions are added to the green wells. Declining green intensity indicates declining NM concentrations (from row 11 to row 2).

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		9/11

Note: Due to the 1:2 dilution of the NMs final concentrations are halved.

- Incubate plate in a humidified incubator at 37°C and 5% CO₂.
- Measure fluorescence in a multi-well plate reader after 1, 2, 3 and 4 hours. After each measurement place plate back into incubator.

Fluorescence settings: excitation at $\lambda=485-495$ nm, emission at $\lambda=517-528$ nm.

6.6 Data evaluation

Data are presented as fluorescence values and represent the mean of three technical replicates and their standard deviation.

Wells B2-D11 (which do not receive DCF dye) serve as control wells to determine NM intrinsic fluorescence. If these values show any concentration dependent trend, they should be subtracted from DCF-containing corresponding values (blank correction as described in SOP “Detection of reactive oxygen species in A549 or THP-1 cells”).

To better illustrate the quenching effect of a given NM, solvent control values (E2-G2) can be normalized to 100%. All other values are then expressed in % of the solvent control and can be interpreted as “quenching efficiency in %”.

7 Quality Control, Quality Assurance, Acceptance Criteria

8 Health and Safety Warnings, Cautions and Waste Treatment

9 Abbreviations

DA	diacetate
DCF	2', 7'-Dichlorofluorescein
ddH ₂ O	double-distilled water
H ₂ DCF	2', 7'-Dichlorodihydrofluorescein
H ₂ DCF-DA	2', 7'-Dichlorodihydrofluorescein-diacetate
HBSS	Hank's balanced salt solution
NM	nanomaterial
ppm	parts per million
RT	room temperature

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		10/11

10 References

Belyanskaya L, Manser P, Spohn P, Bruinink A, Wick P; 2007: The reliability and limits of the MTT reduction assay for carbon nanotubes-cell interaction. *Carbon* 45: 2643-2648

Casey A, Herzog E, Davoren M, Lyng FM, Byrne HJ, Chambers G; 2007: Spectroscopic analysis confirms the interactions between single walled carbon nanotubes and various dyes commonly used to assess cytotoxicity. *Carbon* 45: 1425-1432

Guo L, von dem Bussche A, Buechner M, Yan A, Kane AB, Hurt RH; 2008: Adsorption of essential micronutrients by carbon nanotubes and the implications for nanotoxicity testing. *Small* 4: 721-727

Kroll A, Pillukat MH, Hahn D, Schnekenburger J; 2009: Current *in vitro* methods in nanoparticle risk assessment: Limitations and challenges. *Eur. J. Pharm. Biopharm.* 72(2): 370-377

Monteiro-Riviere NA, Inman AO; 2006: Challenges for assessing carbon nanomaterial toxicity to the skin. *Carbon*44: 1070-1078

Pulskamp K, Diabate S, Krug HF; 2007: Carbon nanotubes show no sign of acute toxicity but induce intracellular reactive oxygen species in dependence on contaminants. *Toxicol. Lett.* 168: 58-74

Wörle-Knirsch JM, Pulskamp K, Krug HF; 2006: Oops they did it again! Carbon nanotubes hoax scientists in viability assays. *Nano Lett.* 6: 1261-1268

Document Type	Document ID	Version	Status	Page
SOP	O_DCF_interference	1.1		11/11